

Benchtop Micro-ESR™



Active Spectrum manufactures the world's smallest Benchtop Electron Spin Resonance (EPR/ESR) spectrometer.

Industrial, scientific and educational applications of electron spin resonance are now extraordinarily cost-effective.

ASI's Benchtop Micro-ESR™ includes our patented miniature electron spin resonance spectrometer, operating at 9.5 GHz with a sweep range of 300 Gauss and sub-micromolar sensitivity.

Also included are an automatic temperature controller and full Windows 7 computer system with Ethernet and multiple USB ports.

Active Spectrum's Benchtop Micro-ESR™ spectrometer is used in industrial and academic labs and educational institutions worldwide. Some sample applications are:

- Crude oil analysis: asphaltene and vanadium content
- Lubricants analysis: oxidation of engine oil, hydraulic oil and turbine oil
- Spin-trapping (PBN, TEMPOL, BMPO, DMPO)
- Shelf life of food products (vegetable oil, beer, wine)
- Biodiesel oxidative stability
- Detection of raw fuel dilution in marine engines
- Additives
- Soot
- Catalyst coking



SPECIFICATIONS

| Parameter | X-band (9.5 GHz) |
|---------------------------|--|
| Sensitivity | 0.1 μ M |
| Sample Tube diameter | ϕ 5.8 mm |
| Magnet Assembly | 3480 \pm 150 Gauss |
| Magnetic Field Uniformity | 0.25 Gauss |
| Supply Voltage | 15 VDC / 6.7A (120V/240V Wall Adapter Included) |
| Data Interfaces | Ethernet and USB |
| Screen | 8.4" Touchpanel display with Windows 7 Embedded DVI/HDMV/VGA |
| Dimensions | 12" x 12" x 12" Case |

SAMPLE RESULTS

Sensitivity Demonstration: Micro-ESR Spectrometer S/N 21.6 with 1 μ M TEMPOL in water

| Spectrometer Settings | Sample |
|---|--|
| Analyte: 1 micro-molar TEMPOL in deionized water. | Sample Cell: Glass flat cell (5mm x 0.5mm) |
| Microwave power: 25 mW | Sample Volume: 50 μ L |
| Sweep Width: 60 Gauss | |
| Sweep Duration: 3 minutes | |
| Modulation field: 1 Gauss | |

Signal-to-Noise Calculation

| Parameter | Value |
|------------------|----------|
| S/N (Amplitude) | 21.60 |
| Signal Amplitude | 21.17 mV |
| Noise Amplitude | 0.98 mV |
| S/N (Power) | 466.6 |




Figure 1: 3-minute scan of 1 micro-molar TEMPOL in aqueous solution

FEATURES



Data Interface



Sample Insertion



Touchpanel Display

AVAILABILITY

Please contact us at +1 650-212-2625 or sales@activespectrum.com to discuss your application in more detail.