



Analyzing Spin-labeled ApoA-I via SDSL with Micro-ESR

Elizabeth E. McCarthy^{1,2}, Colin Elliott², and Dr. James R. White²
(1) University of California, Berkeley, Department of Chemistry, Berkeley, CA, 94720
(2) Active Spectrum, Inc., Foster City, CA, 94404



54th Rocky Mountain Conference on Analytical Chemistry, July 15th - 19th, 2012.

1 Abstract

Background

- High density lipoprotein (HDL), "good" cholesterol, helps prevent coronary artery disease (CAD).
- Micro-ESR can measure functional HDL levels via site directed spin labeling of the protein apolipoproteinA-I (apoA-I).

Micro-ESR

- Sensitivity
 - Detection of spectral differences that depend on spin label location on the protein.
- Quantitative Ability
 - Use of Manganese reference to provide quantitative data.

2 Background

Atherosclerosis is a major cause of coronary artery disease (CAD). High density lipoprotein (HDL), commonly known as the "good" cholesterol, serves a protective role through preventing atherosclerosis by removing cholesterol from the body. Thus, low functional HDL levels can raise the risk of CAD. ApoA-I is an important protein in high density lipoprotein (HDL). Thus, if the functional HDL levels are abnormally low, the risk of atherosclerosis and thus CAD increases.

3 SDSL

Site Directed Spin Labeling (SDSL)

- Uses electron spin resonance spectroscopy (ESR) to determine protein structure
- Uses a stable nitroxide as a spin label on the target residue
- Line shape and width of the peaks of the ESR signal reveal the structure of the protein near the residue

Broader Signal 

- Ordered protein structure which allows for less rotation of the residue
- Dipolar coupling

Narrower Signal 

- Less structured protein which allows for more rotation of the residue

4 Comparison of ESR Spectra of Spin-Labeled ApoA-I




Figure 1. Micro-ESR spectrum plots of spin labeled ApoA-I. The differences in the rotational freedom and dipolar coupling of the spin label at the different sites are reflected in the line shape and width of the spectrum plots.

References

- Bordignon, Enrica. 2010. Nitroxide Spectrum Analysis Tutorial. http://www.epr.ethz.ch/news/Bordignon_tutorial_Nitroxide_spectrum_analysis.pdf. Web. 17 Apr. 2012
- Lagerstedt, Jens O., Madhu S. Budamagunta, Michael N. Oda and John C. Voss. 2007. "EPR Spectroscopy of Site-Directed Spin Labels Reveals the Structural Heterogeneity in the N-Terminal Domain of ApoA-I in Solution." *Journal of Biological Chemistry*. 282:9143-49.
- Oda, Michael, Trudy Forte, Robert Ryan, and John Voss. 2003. "The C-terminal domain of apolipoprotein A-I contains a lipid-sensitive conformational trigger." *Nature Structural Biology*. 10(6):455-460.

6 ESR Spectrum for Quantitative Analysis

G217C with MnO reference




Figure 2. Micro-ESR spectrum plot of G217C residue on apoA-I with manganese oxide as an external reference. ApoA-I signal can be quantified via comparison to the reference signal (the outermost peaks on either side).




Figure 3. Micro-ESR Spectrometer (12"x12"x12", 13.8 lbs.)

7 Acknowledgements

- Colin Elliott, Active Spectrum, Inc.
- Dr. James White, Active Spectrum, Inc.