

Rectangular ESR Sample Tubes

- Get more than double the signal to noise on aqueous samples with rectangular sample tubes
- Active Spectrum's rectangular sample tube holder keeps tubes perfectly aligned
- Tubes are easily filled with a pipetor with gel tips
- Will work with any polar solvent



Available from Active Spectrum:

- | | |
|---|----------------|
| • Rectangular Sample Tubes (set of 5) | \$50.00 |
| • Rectangular Sample Tube Holder | \$300.00 |
| • Rectangular Tube Starter Kit | \$400.00 |
| • Extended Rectangular Tube Accessory Kit | Call for quote |



240:1 -1.7 mm capillary tube

580:1 -rectangular tube



Both tubes were filled with the same aqueous TEMPOL sample, and run under identical conditions:

16 scans, 15 mW