


# NEW


## Micro-ESR with Extended Range



- 100-4000G minimum sweep range
- Compact bench top design
- Competitive pricing
- Variable temperature accessories available



VO(acac)<sub>2</sub> in CHCl<sub>3</sub>



PBN-DMSO-OH spin adduct

# Technical Data

<b>Frequency</b>	9.6 to 9.8 GHz
<b>Resonator</b>	Cylindrical Dielectric
<b>Sample Tube Diameter</b>	Up to 5.8 mm
<b>Sweep Range</b>	Up to 6500 Gauss (water cooled); Up to 4000 Gauss (air cooled)
<b>Supply Voltage</b>	120/240 VAC
<b>Data Interfaces</b>	Ethernet, USB, and 802.11 WiFi
<b>Screen</b>	21 cm Touch Panel display with Windows 7 Embedded DVI/HDMI/VGA
<b>Dimensions</b>	30.5 x 35.5 x56 cm <sup>3</sup>
<b>Mass</b>	54 kg